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1. Intended purpose of the product

The elastic tape is intended for application to intact skin in order to provide functional support for muscles and joints by stabilizing soft tissues, supporting proprioception, and aiding the body's natural regenerative processes.

2. Indications for use

- The product may be used:
- to support muscle and joint function
 - to stabilize soft tissues
 - to improve proprioceptive sensation
 - as an aid during periods of musculoskeletal overload
 - as an aid during physical activity and throughout the recovery process

3. Contraindications

- Do not use:
- on damaged, irritated, or diseased skin
 - on wounds, abrasions, burns, or active skin lesions
 - in case of hypersensitivity or allergy to the adhesive ingredients or the tape material
 - if skin reactions have occurred after previous use of the product

- in cases of skin diseases such as eczema, psoriasis, or other chronic dermatological conditions
- if skin reactions have occurred after previous use of the product
- in cases of skin diseases such as eczema, psoriasis, or other chronic dermatological conditions
- in cases of venous thrombosis or other circulatory disorders without prior consultation with a doctor
- in cases of fractures, including improperly healed fractures, recent injuries, acute inflammatory conditions, or fresh postoperative scars without prior medical consultation
- directly after getting a tattoo
- in cases of varicose veins
- in cases of diabetes
- in cases of aneurysms

Use in children

- the tape may be used in children only under adult supervision and after consultation with a doctor or physiotherapist
- particular attention should be paid to the sensitivity of a child's skin, and prolonged use should be avoided

If there is any doubt regarding the use of the product, consult a doctor or physiotherapist.

4. Warnings and precautions

- the product is intended for external use only
- apply only to clean, dry, and intact skin
- do not apply the tape with excessive tension
- if itching, burning, redness, pain, or any other signs of irritation occur, remove the tape immediately
- do not apply in a way that impairs proper blood circulation.

- do not reuse the product after removal
- the product should be used in accordance with a specialist's recommendations or as intended
- keep the product out of the reach of children
- all serious incidents involving the product must be reported to the manufacturer and to the competent authorities of the Member State in which the user or patient is located.

5. Directions for use

1. make sure that the skin at the application site is clean, dry, and free from oil
2. if necessary, cut the tape to the appropriate length
3. round the corners of the tape to improve adhesion
4. remove the protective backing paper and apply the tape to the skin according to the intended method of application
5. apply the ends of the tape without tension
6. after application, gently rub the surface of the tape to activate the adhesive layer
7. the wearing time should be adjusted individually depending on the condition of the skin and the intended purpose of use

6. Removal of the product

Remove the tape slowly and carefully, in the direction of hair growth.

For sensitive skin, it is recommended to remove the tape after moistening it with warm water or oil.

7. Storage

Store in a dry place at room temperature, away from sources of heat, moisture, and direct sunlight.

8. Package contents

1 x elastic tape

9. Final notes

The product does not replace medical consultation or specialist treatment.

If symptoms persist or worsen, consult a doctor.

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