**Instruction Manual and Warranty for the PRO THERAPY
MAX Variable Amplitude Massage Gun**

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
PL POLSKI

Strona 51

Thank you for purchasing our muscle massager with adjustable amplitude.

Before using this product, please read this manual carefully and keep it for future reference.

Read this manual carefully before use and pay special attention to the safety warnings.

Explanation of symbols			
	Direct current		Li-ion recyclable
	WEEE Directive identification		Do not roll
	WARNING! Refer to Documents		Please read the manual carefully
	Fragile		Upwards
	Keep dry/ Protect from rain		Paper recycling
	Stacking layer limit		

Safety Instructions

Before using the device, please read and understand the following safety instructions to avoid injury or danger.

About the User Groups and Body Parts

- For use only by healthy adults.
- Use only on soft tissue areas of the body, avoiding pain or discomfort.
- Do not use on the head or on any hard or bony parts of the body.
- Bruising may occur even under controlled use. Check the treated area frequently and stop using the device immediately if pain or discomfort occurs. If symptoms persist, seek medical attention.
- Use only on clean, dry skin. Apply gentle pressure and move the device slowly over the skin. Do not keep it in one area for an extended period.
- Keep fingers, hair, and other body parts away from the moving parts or back of the massage head to prevent pinching or injury.
- Avoid using the device on areas with little muscle tissue or near bones, such as the spine, neck, joints, tibia, back of the hand, or top of the foot.
- Do not use near sensitive areas, organs, or the central nervous system, including the eyes, teeth, breast implants, prostheses, chest, abdomen, spine, or reproductive organs.
- Do not use on or near areas containing implants, as this may affect their function or cause loosening.
- Do not use on areas with bleeding, unhealed wounds, local infections, or inflammation.
- People with conditions that reduce sensitivity or the ability to feel pressure should not use this product.
- People with circulatory or thrombotic disorders should not use this product.
- People who are pregnant or have recently undergone surgery should not use this product.
- People with heart failure or severe weight loss (cachexia) should not use this product.

Important Safeguards

- Do not use the product on any area where a bone fracture may be present.
- If pain persists, stop using the device immediately and consult a doctor.
- This product is not a medical device and is not intended for medical use.
- Check the device and its connecting cable before each use. If any damage is found, it must be replaced by the manufacturer or an authorized service provider to prevent injury or danger.

Important

This product should not be used under any of the following conditions without prior approval from a doctor: Pregnancy; complications of diabetes (such as neuropathy or retinal damage); presence of a pacemaker; recent surgery; epilepsy or migraine; osteoporosis; herniated disc; spondylolisthesis; spinal joint stiffness; recent joint replacement; intrauterine device (IUD); presence of metal implants such as needles, plates, or pins; or any other medical condition affecting your physical health.

These contraindications do not necessarily mean that you cannot use this product; however, you must consult your doctor and obtain professional medical advice before use.

About Device Operation

- Select the appropriate massage head before using the device.
- Do not place any objects into the ventilation openings.
- Do not immerse the device in water or allow liquid to enter the ventilation holes.
- Do not drop or misuse the device.
- Check the device, battery, and charger before each use.

Important Safeguards

- Do not modify or use the device for any other purpose.
- Do not leave the device running or charging unattended.
- Do not operate the device while it is charging.
- Do not throw the device into fire.
- Do not operate the device continuously for more than one hour. Allow it to rest for at least 30 minutes before the next use.

About battery

- Do not throw the battery into water or allow it to get wet.
- Do not charge the battery in hot environments or near open flames or heaters.
- If you notice battery leakage or a strong odor, move the device away from heat or flames immediately.
- Do not throw the battery into fire or expose it to high temperatures.
- Do not transport or store the battery with metal objects such as necklaces, hairpins, or coins.
- Do not connect the positive and negative terminals of the battery with wires or metal objects to avoid a short circuit.
- Do not pierce the battery with nails or other sharp objects.
- Do not strike, step on, crush, or throw the battery.
- Do not disassemble, modify, or attempt to open the battery in any way.
- Do not place the battery in a microwave oven or pressure cooker.
- If you notice a strong odor, excessive heat, deformation, discoloration, or any other abnormal condition, stop using the battery immediately.
- If the battery is in use or charging, disconnect it from the device or charger immediately.

Product Overview


This muscle massager is designed to provide a comfortable and relaxing massage experience. It helps ease muscle tension, supports relaxation after physical activity, and promotes a general feeling of comfort and well-being.

Purpose


For physical massage after exercise.

Main Technical Specifications

- Product name: PRO THERAPY MAX variable amplitude muscle massager;
- Product model: T-POZBEO001;
- Structure: The main structural components include the control unit, massage heads, and lithium battery;
- Massage gun vertical vibration frequency: ≤ 60 Hz;
- Amplitude: 8-16mm;
- Device: Rated voltage: 25W;
- Instrument: rated voltage: 21.9V $\approx$;
- Battery capacity: 2600mAh;
- Adapter: Input: 100-240V~50/60Hz 0.5A
Output: 5V $\approx$ 3A/9V $\approx$ 2.22A/12V $\approx$ 1.67A
- Intelligent auto-off: 10 minutes;
- Speed levels: Gear 1:1500r/min; Gear 2:1700r/min;
Gear 3:1900r/min; Gear 4:2100r/min;
Gear 5:2300r/min; Gear 6:2500r/min;
- The device is equipped with a switch to control how it works;
- The massage head can be replaced as needed.

Main Structure

Main Structure


Installing and Replacing the Massage Head

Before replacing a massage head, make sure the device is switched off and unplug the massage head from the installation slot. Gently apply force to detach the massage head. Attach the new head to the massage gun by inserting the tab into the slot and ensure it is secured.

 Air Cushion Head Ideal for relaxing and massaging sensitive muscles of the body.	 Round Head Ideal for all muscle groups.
 Bullet head Ideal for targeting smaller muscle groups and precise deep-tissue relief.	 Flat head Suitable for all muscle groups. Provides an even and comfortable massage.
 Spade head Ideal for shoulders blades, IT bands, scraping, and flushing techniques.	 Dome Head Ideal for both small and large muscle groups, offering versatile full-body use.
 Ball Head Ideal for large muscle groups such as the quadriceps and gluteus maximus. Ideal for relaxing tight muscles and relieving general soreness.	 U-Shaped Head Ideal for the spine, the neck, shoulders, and other sensitive muscles and soft tissues.

Operation Steps

- 1) Make sure the device is turned off before attaching the massage head.
- 2) Press the massage head firmly into the device. (Tip: For the U-shaped head, make sure the buckle is aligned with the slot before inserting it. Otherwise, the head may fall off during use.)

3) Button function description refers to the following table.


Button	Operation	Function
⏻	Long press	Turn on the device while it is off
		Turn off the device while it is on
	Press	In M1 mode or standby state, increase the gear In the standby state of M2/M3/M4 mode, the device works
⏮	Long press	In standby mode, mode switching In the working state of M2/M3/M4 mode, enter the standby state
		In the working state of M1 mode, reduce the gear
	Press	In the working state of M2/M3/M4 mode, enter the standby state
⏪	Rotation	"+" amplitude when the device is in standby or running
⏩	Rotation	"-" amplitude when the device is in standby or running

Note:

1. When the device is powered on, it defaults to M1 mode.
2. Mode switching can only be performed in standby mode.
3. M1 is a common gear mode (0-6 levels adjustable), while M2, M3, M4 are fixed modes and cannot be adjusted.
4. The device will automatically enter the standby state after operating for 10 minutes.
5. The device automatically turns off after 60 seconds of inactivity in standby mode.


Operation

4) Display Description

5) The arm angle of the device can be adjusted within a range of 0-90°.

To rotate the arm, press and hold the arm angle adjustment switch located on the left side of the device.

Operation

Mode Explanation

Mode	Introduction	Application	Features
M2	Automatic incremental mode. 1500, 1700, 1900, 2100, 2300, and 2500 rpm work for 5 seconds respectively, cycle and repeat. The OLED screen will show the speed it's in. The massage gun will automatically turn off in about 10 minutes.	Muscle relax and recovery after workout	Low initial speed, speed will change from lowest to highest and repeat
M3	Automatic decreasing mode. 2500, 2300, 2100, 1900, 1700, and 1500 rpm work for 5 seconds respectively, cycle and repeat. The OLED screen will show the speed it's in. The massage gun will automatically turn off in about 10 minutes.	Activate muscle before workout	High initial speed, speed will change from highest to lowest and repeat
M4	Triangle wave mode. 1700, 1800, 1900, 2000, 2100, 2200, 2300, 2200, 2100, 2000, 1900, 1800 work for 2 seconds respectively, cycle and repeat. The OLED screen will show the speed it's in. The massage gun will automatically turn off in about 10 minutes.	For all needs	The speed will change rapidly

Recommended usage method

THIN MUSCLES SHALLOW MASSAGE		
 Shoulder muscle group Amplitude: 8-9mm 1500-1900r/min Gear 1-3	 Chest muscle group Amplitude: 8-9mm 1500-1900r/min Gear 1-3	 Back muscle group Amplitude: 8-9mm 1500-1900r/min Gear 1-3
THICK MUSCLES DEEP MASSAGE		
 Calf muscle group Amplitude: 10-12mm 1900-2100r/min Gear 3-4	 Leg muscle group Amplitude: 10-14mm 1900-2300r/min Gear 3-5	 Gluteal muscle group Amplitude: 14-16mm 2100-2500r/min Gear 4-6

Charging

- We recommend using the power adapter provided with the product to ensure optimal and stable charging. The device may also charge with other compatible adapters, but charging time and performance can vary. If the “” icon appears on the screen, please use the supplied adapter or replace the adapter being used.
- Charge the device for approximately 3 hours before first use.
- Connect the output end of the charging cable to the TYPE-C charging port located on the bottom of the product. Then connect the input end of the charging cable to the power adapter.
- When the “” icon appears on the screen, and the charging indicator light is green and flashing, the device is charging.
- When the “” icon appears and the charging indicator light is off, the device is fully charged.
- When the “” icon appears on the screen, it means the battery level is low and the device needs to be charged.
- When the “” icon appears on the screen, it indicates battery overheat protection.
- Please let the device rest for at least 30 minutes before using it again.
- When the “” icon appears on the screen, it means a 5V power adapter is being used.
- When the “” icon appears on the screen, it means a 9V power adapter is being used.
- When the “” icon appears on the screen, it means a 12V power adapter is being used.
- It is not recommended to fully discharge the battery before charging.

Troubleshooting

Problem	Cause	Solution
The massage gun cannot be turned on.	The power is off.	Press and hold the Power button to turn on the massage gun.
	Low battery.	Charge the massage gun.
	The massage gun automatically shuts down after 10 minutes.	Press and hold the Power button to continue using the massage gun for another 10 minutes.
The massage gun automatically shuts down.	Overload protection is activated; the massage gun will shut down automatically when overloaded.	Turn on the massage gun again after 10 mins.
	The massage gun automatically shuts down after 10 minutes	Turn on the massage gun by pressing and holding the Power button to use it for another 10 minutes.
The massage gun is blocked.	Lubricating oil may solidify due to low temperature.	Loosen the massage head slightly, then turn on the massage gun.

Maintenance

- Clean the device with a slightly damp cloth, then wipe it dry.
- Do not allow any liquids to enter the device.
- Store the device in its box when not in use for an extended period.
- Do not disassemble the device.
- To ensure optimal battery performance, fully charge it to 100% and discharge it to 0% three times after purchase.
- If the device is not used for a long time, recharge the battery once per month.
- Damage caused by unauthorized disassembly or repair is not covered by the warranty.

Waste Battery Disposal

Please support the recycling of rechargeable batteries (lithium-ion batteries).

 Danger	The battery used in this product is a special rechargeable battery and must not be used in other devices. Do not remove the rechargeable battery for charging. Failure to follow these instructions may cause overheating, fire, explosion, or damage. • Do not throw the battery into fire or expose it to high temperatures. • Do not pierce, strike, crush, or deform the battery. • Do not disassemble or modify the battery. • Do not connect the positive and negative terminals of the battery with metal objects. • Do not store the battery with metal items such as necklaces, hairpins, or coins. • Do not charge, use, or store the battery near heat sources or in direct sunlight. Failure to observe these precautions may result in overheating, fire, explosion, or other serious hazards.
 Danger	• Keep removed rechargeable batteries out of the reach of children. • If a battery is accidentally swallowed, seek medical attention immediately. • If battery fluid leaks, avoid contact with skin and handle as follows: – If the fluid gets into the eyes, it may cause blindness. Do not rub your eyes. Rinse immediately with clean water and seek medical attention. – If the fluid contacts the skin or clothing, it may cause irritation or burns. Wash the affected area with clean water. If discomfort persists, seek medical advice. • Disassembly is only permitted when the product is being disposed of. Failure to follow this instruction may result in fire, electric shock, or injury. • For repairs (other than battery replacement), contact the original retailer or an authorized service center.

- The built-in lithium-ion battery is a valuable resource that can be recycled.
- Please have it removed and taken to a qualified recycling or waste-disposal facility.

Conditions for Storage

- Store the product in a dry, well-ventilated, non-corrosive environment.
- If the device will not be used for an extended period, clean it, pack it properly, and store it.
- Recharge the device once a month to maintain battery performance and prevent moisture damage.

Warranty Conditions

1. For the purposes of the guarantee:
 - a) "Manufacturer" refers to HABYS Spółka z ograniczoną odpowiedzialnością (Limited Liability Company) with its registered seat in Jasło, ul. Produkcyjna 16, 38-200 Jasło, entered into the Register of Entrepreneurs kept by the District Court in Rzeszów, XII Commercial Division of the National Court Register, with the assigned KRS number: 0000513317, with the share capital in the amount of PLN 10.750.000 (ten million seven hundred fifty thousand zlotys), the capital has been fully paid up, NIP: 6852208438, REGON: 180186290.
 - b) "Buyer" refers to the final Buyer of the equipment, manufactured by the Manufacturer, including both the consumer and the entrepreneur.
 - c) "Distributor" refers to an entity that, under an agreement with the Manufacturer, sells the Equipment in a designated area.
 - d) "Authorized Service Centre" refers only to a service held by the Manufacturer or a service indicated by the Manufacturer and held by its Distributor.
 - e) "Equipment" refers to a movable item, manufactured by the Manufacturer, being the subject of the sales and to which the said warranty applies.

2. The manufacturer ensures good quality and efficient operation of the equipment for which the warranty card has been issued, for the period of 2 years.
3. The guarantee repair does not include the activities provided for in the manual, which the Buyer is obliged to perform on his own (installation of the Equipment, maintenance of the Equipment, etc.).
4. Liability under the guarantee covers only defects resulting from causes inherent in the Equipment., In particular, the guarantee does not cover any damage resulting from improper or inconsistent use, storage, maintenance or transport of the Equipment, as well as colour differences between the parts exposed and not exposed to solar radiation, the influence of chemical or abrasive agents, atmospheric pollution or accidental punctures and abrasions, the presence of stains or rings from greasy or dyeing substances (fats, creams, lotions or oils), permanent staining with dyes such as inks, permanent markers or other non-permanent dyes used in the production of clothing (e.g. dyes used to make blue jeans).
5. The guarantee does not cover any mechanical damage to the Equipment.
6. The buyer loses the guarantee rights in the event of:
 - a) the Manufacturer states that the Equipment was

- repaired by entity other than the Authorized Service,
 - b)making any design modifications to the Equipment without the prior written consent of the Manufacturer,
 - c)using the Equipment contrary to its intended use,
 - d)using the Equipment after revealing any defect in it.
- 7.Regardless of the loss of the guarantee rights, described above, in point 6, the Manufacturer shall not be liable for any damage to property or caused to third parties, which may be caused by the use of the Equipment covered by the guarantee.
- 8.The guarantee does not exclude or limit the rights resulting from the warranty for defects. Exercising the rights under the guarantee does not affect the seller's liability under the warranty for defects. If the Buyer exercises the rights under the guarantee, the time limit for the exercise of the rights under the warranty for defects is suspended on the date the Manufacturer is notified of the defect. The period runs further from the date of the Manufacturer's refusal to perform the obligations arising from the guarantee or the ineffective expiry of the time for their performance.
- 9.The scope of the guarantee coverage is limited to the territory of the country in which the Equipment

- was sold.
- 10.The Manufacturer shall perform its obligations under the guarantee (removal of defects or damage to the structure of the Equipment arising from reasons inherent in the Equipment, and disclosed during the guarantee period) free of charge within 14 days from the date of delivery of the Equipment to the Manufacturer or the Authorized Service Centre. If the repair of the Equipment cannot be performed within the above-mentioned period, the Manufacturer shall inform the Buyer of the above, and also indicate the reasons for the delay and the expected date of fulfilment of the recognized guarantee request.
- 11.The use of the guarantee rights requires the total (cumulative) fulfilment of four conditions:
- a)presenting a VAT invoice, receipt or other proof of purchase of the equipment to the Manufacturer, photographic documentation of the defective Equipment, together with a brief description of the reasons for submitting a guarantee claim,
 - b)reporting the defect within 14 (fourteen) days of its detection,
 - c)delivering the Equipment in original packaging at the expense of the Manufacturer, to the address indicated above in the guarantee card or to the

- address of the Authorized Service Centre,
d) designation of the place where the Manufacturer delivers the equipment after repairing it (limited to the country in which the Equipment was purchased).
12. The guarantee period is extended for the duration of the repair, starting from the date of delivery of the Equipment to the Manufacturer or to the Authorized Service Centre, to the date on which the Equipment was returned to the Buyer.
13. The responsibility of the Manufacturer or the Authorized Service Centre is to define the scope and the method of repairing the Equipment.
14. The guarantee is limited to the repair of the Equipment. The Manufacturer reserves the right to replace the Equipment with a non- defective one, if the repair exceeds the price of the Equipment indicated in the proof of purchase. Under the guarantee, the Buyer shall not be entitled to demand that the Equipment be replaced with a new one. The guarantee, in particular, does not cover liability for losses or damages incurred by the Buyer in connection with the fact that one is not able to use the Equipment.
15. The Manufacturer shall not be liable for any damage to the Equipment caused during the performance of the transport service by the delivery

- company until it is collected by the Manufacturer from the delivery company.
16. If the shipment of the Equipment was ordered by the Manufacturer, the Buyer is obliged to check the technical condition of the Equipment upon its receipt. If the Buyer finds that the Equipment has been damaged during the transport, the Buyer is obliged to:
- a) immediately, not later than within 3 days from the date of receipt of the shipment, notify the Manufacturer of any damage to the Equipment that occurred during transport, in order to enable the Manufacturer to submit a complaint to the delivery company and pursue related claims,
- b) provide the Producer with a complaint report drawn up between the Buyer and the delivery company.
17. The guarantee regulations contained in this document are the sole and exclusive guarantee rules for HABYS limited liability products.

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